INDIRECT FLUORESCENCE ASSAY FOR HUMAN HERPESVIRUS 6 (HHV6) IgM ANTIBODY

EXPORT ONLY

I-HV601M	120 Tests
I-HV602M	40 Tests
I-HV603M	40 Tests

Intended Use

The SCIMEDX Indirect Fluorescence Assay (IFA) for Human Herpesvirus 6 (HHV6) IgM Antibody is intended for the qualitative and semi-quantitative detection of IgM (Immunoglobulin M) antibody to HHV6 in human serum or plasma. Detection of HHV6 IgM antibody in humans can be used as an aid in the diagnosis of primary infection with this virus.

Introduction and Summary of Test Procedures

Human herpesvirus 6 (HHV6) is a human lymphotropic virus, first isolated from immunosuppressed individuals and patients with lymphoproliferative disorders. It is a DNA virus, originally named human B-lymphotropic virus (HBLV), now classified as a member of the family of human herpesviruses (1-2). HHV6 is genetically and serologically distinct from the other known human herpesviruses (2-4). Serological studies, using the indirect fluorescence assay (IFA), have shown that antibody to the virus is common in the general population. Primary infection with HHV6 usually occurs before the age of two. A high percentage of humans above the age of one are seropositive for anti-HHV6 antibody (4-6). Primary infection in adults is rare because of the high rate of infection in childhood. While primary infection in healthy adults usually results in mild illness, more serious complications can arise in immunocompromised patients or patients undergoing organ or bone marrow transplantation (4). HHV6 has been identified as the causative agent of roseola infantum, also known as exanthem subitum. The virus has been isolated from children with exanthem subitum and the association was further supported by serologic studies (7-9). Roseola infantum is considered a benign non-seasonal disease of childhood with clinical indications of rash and fever. Infants are protected by maternal antibodies so that infection is most frequent between 6 and 18 months of age. Serological information is valuable in diagnosis since detection of antibody by the IFA technique is both sensitive and specific, without cross-reaction to other human herpesviruses (3, 5, 10, 11-13).

Principle of the Test

SCIMEDX fluorescent antibody assays use the indirect method of antibody detection and titer determination. Patient serum or plasma samples are applied to cultured cells containing inactivated viral antigens provided on paint delineated wells on glass microscope slides. During a 60 minute incubation, antibody specific for HHV6 antigens forms an antigen/antibody complex with the HHV6 antigens in the infected cells. In a brief washing step, nonspecific antibodies and other unreacted serum proteins are eliminated. Fluorescein-conjugated goat anti human IgM is then applied to the wells of the glass slide. The anti-IgM conjugate combines with human IgM, if present, during a 30 minute incubation. After a brief wash to remove unreacted conjugate, the slides are viewed by fluorescence microscopy. A positive antibody reaction is denoted by bright green fluorescence at the antigen sites.

Materials Furnished and Storage Conditions

HHV6 Antigen Slides: Slides of human lymphocytes infected with HHV6 on each glass well. The slides are ready for use after removal from protective pouch. Store at 2-8°C. Slides are stable at this storage condition until the expiration date stated on the pouch label.

HHV6 IgM Positive Control: Each vial contains 0.5 ml HHV6 IgM antibody positive human control. This component is a ready for use liquid. Store at 2-8°C. Liquid positive control is stable at this storage condition until the expiration date stated on the vial label.

HHV6 IgM Negative Control: Each vial contains 0.5 ml HHV6 IgM antibody negative human control. This component is a ready for use liquid. Store at 2-8°C. Liquid negative control is stable at this storage condition until the expiration date stated on the vial label.

Fluorescein Conjugate: Each vial contains 1.5 ml fluorescein conjugated goat (inactivated) antihuman IgM (μ chain specific) with Evans Blue and Rhodamine counterstains. The fluorescein conjugate is a conjugation of affinity chromatography purified anti-human IgM with fluorescein isothiocyanate (FITC). Adding Evans Blue and Rhodamine counterstains to the conjugate masks nonspecific fluorescence of the tissue culture cells. This component is a ready for use liquid. Store at 2-8°C. Liquid conjugate is stable at this storage condition until the expiration date stated on the vial label.

Coverslip Mounting Media: Each vial contains 2.0 ml phosphate buffered glycerol with fade retardant. This component is a ready for use liquid. Storage temperature may range from refrigerated to room temperatures (2-30°C). Mounting media is stable at either storage condition until the expiration date stated on the vial label.

Phosphate Buffered Saline (PBS): Each aluminized sealed packet of powdered buffer makes one liter of 1X PBS. Storage temperatures may range from refrigerated to room temperature (2-30°C). Add the entire contents of a PBS packet to one liter of freshly prepared distilled or deionized water. Note: Addition of the salts while rapidly stirring the water will facilitate solubilization. Store PBS as a solution at 2-8°C.

Special Blotters: Absorbent blotters have pre-cut holes for use in drying the slide mask. Storage temperatures may range from refrigerated to room temperature (2-30°C).

General Precautions

IFA Test Kit: No US Standard of Potency.

- Store all kit components at their recommended or suggested temperatures. **Do not freeze.**
- Do not use the components beyond the stated expiration date of each component.
- SCIMEDX optimizes all of the active components in each lot of its IFA kits as a unit. Do not mix components from different lots or from different sources.
- The controls and conjugate contain 0.095% sodium azide that, if allowed to accumulate, can form explosive compounds in lead and/or copper plumbing. Flush drains thoroughly if used to dispose of these materials.
- Use separate pipette tips for each sample and reagent to avoid cross contamination.
- Reagents should be inspected for evidence of bacterial or fungal contamination.

Antigen slides: All IFA antigen slides have fixed cells that do not contain any viable infectious agents. However, good laboratory practices (GLP) require careful slide handling and disposal as with any other potentially biohazardous laboratory material.

- Do not remove the slides from their protective pouch until ready for use.
- Do not reuse substrate slide.

Human Controls: The human controls in these kits have all been tested for hepatitis B surface antigen (HBsAg) and human immunodeficiency virus (HIV) antibody by FDA licensed methods and found to be non-reactive. However, no test system can ensure the absence of these agents. Handle all human serum components, including those received in your laboratory for testing, as potentially biohazardous.

Xn-Harmful Substance

Control and Conjugate Safety Precaution:

Concentration of sodium azide in these components is classified as Harmful and subject to the following risk phrase: "Harmful if swallowed and contact with acids liberates very toxic gas."

US

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention: • Wash thoroughly after handling. • Wear eye/face protection Response: • IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. • If eye irritation persists: Get medical advice/attention. Storage/Disposal: • Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.
Pictogram		
Signal Word	WARNING	
Hazard Statement	Causes serious eye irritation.	

EU

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention: • P264 Wash thoroughly after handling. • P280 Wear protective gloves and clothing. Response: • P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing. • P337+P313 If irritation persists, get medical advice/attention.
Pictogram		
Signal Word	WARNING	
Hazard Statement	H319 Causes serious eye irritation	

Specimen Collection, Storage and Limitations of Test Samples

- Separate aseptically collected serum or plasma from the red blood cells and store frozen (-10°C or lower) until ready for testing. Avoid repeated freezing and thawing.
- If desired, store fresh liquid serum or plasma samples at 2° to 8°C for up to one week without loss of antibody activity.
- Do not use excessively lipemic samples without delipidization.
- Do not use contaminated samples.
- SCIMEDX recommends screening dilutions of 1:10 and 1:40 for testing purposes unless a pretreatment step removes interfering IgG. If serum is pretreated, testing at the 1:10 dilution is sufficient. We have available a pretreatment reagent. Please call customer service for additional information.

Additional Materials Required

- Test tubes, racks, pipettes, microtiter plates and safety pipetting devices for making sample dilutions
- 37°C incubator
- Moist chamber for incubating slides
- Slide holder rack and staining dish for washing slides
- Coverslips: 22 X 50 mm No. 1 thickness glass

Fluorescence microscope: A fluorescence microscope equipped with the following was used to calibrate the controls and conjugate:

- 10X eyepiece
- 16X or 40X objectives
- Epi-illuminator with 50W halogen lamp
- FITC-excitation filter KP490
- Yellow absorbing filter K530
- Red suppression filter BG38

The fluorescein label has an excitation peak of 490 nm and an emission peak of 520 nm. Differences in endpoint reactivities and fluorescence intensities may be due to the type and condition of the fluorescence equipment used in your laboratory.

IFA Procedure

1. For IgM antibody determination, prepare a 1:10 and 1:40 screening dilution of each test sample in PBS. Prepare all dilutions in a minimum volume of 0.10 ml with PBS as the diluent. For samples pretreated to remove IgG a 1:10 dilution is sufficient. SCIMEDX can provide a pretreatment reagent to remove IgG by the following procedure: Add 45 μ l of the reagent to 5 μ l of test sample and mix well. The resulting 1:10 dilution can be used immediately. Precipitation may occur but will not interfere with the IFA test results.
 2. Remove slides from protective pouch and apply 1 drop (approximately 20 μ l) of the diluted test sample(s) to each well. Add sufficient volume to completely cover each well, but cross-mixing of contents between wells should not occur.
- Note: Each day's test run requires one well each for positive control, negative control, and PBS (conjugate control).
3. Incubate the slides in a moist chamber for 60 minutes at 37°C.
 4. Rinse the slides in a light stream of buffer. Avoid directing the stream at the wells.
 5. Wash the slides for 10 minutes in PBS, changing the PBS solution after 5 minutes. Agitate the slides by moving the rack up and down in the buffer.
 6. Blot the paint mask surrounding the test wells with the special blotters provided in the kit.
 7. Apply one drop of the ready to use conjugate to each test well.
 8. Incubate the slides in a moist chamber for 30 minutes at 37°C.
 9. Repeat steps 4 (PBS rinse), 5 (10 minute PBS wash), and 6 (blot).
 10. Apply the glycerol mounting media and 22 X 50 mm glass coverslip to cover the wells of the slide.
 11. Observe the reactivity under fluorescence microscopy using 20-40X magnification. For best results, examine slides immediately after completion of the test. To obtain equivalent results, seal slides or keep humidified to minimize dehydration of mounting medium, store in dark at 2-8°C, and read within three days. Positive reactivity may range in fluorescence intensity from brilliant to weak. Grade the fluorescence reaction according to the following intensity scale: 4⁺ (brilliant), 3⁺ (bright), 2⁺ (moderate) and 1⁺ (weak).

Interpretation of Results

- Bright green fluorescent staining of the infected cells denotes an HHV6 IgM antibody positive reaction. The fluorescent staining pattern of HHV6 infected cells is variable. Depending on the cell's stage of infection, the fluorescent pattern can vary from a small portion of the infected cell fluorescing to the whole cell fluorescing. Fluorescence can also range from granular to homogeneous. To provide an internal control, each well on the microscope slide contains both HHV6 infected and uninfected cells. Preparation of the slide in this manner is intentional. Uninfected cells, stained red by the counterstain, provide a contrasting background. Infectivity of the cells ranges from 20% - 60%. Titration of positive HHV6 IgM sera is not necessary except to minimize interfering antibodies or to provide quantitative information. In a titration series, the highest serum dilution demonstrating a 1+ reaction is interpreted as the endpoint.
- Absence of specific fluorescent staining of the infected cells denotes an HHV6 IgM antibody negative reaction.

Significance of Interpretation

1. No discernible fluorescence of the infected cells found at the screening dilutions.	1. Test sample is HHV6 IgM antibody negative.
2. Specific positive fluorescence of the infected cells found at any of the screening dilutions.	2. Test sample is HHV6 IgM antibody positive, indicating current infection.
3. Fluorescence found in both infected and uninfected cells.	3. Test sample is exhibiting a non-specific reaction.

Quality Control

- To ensure the test is working properly, use the positive and negative controls at least once for each day's testing.
- The type and age of the fluorescence microscope and hours of UV bulb usage can affect fluorescence intensity and titration endpoints to some degree. The HHV6 antibody positive control furnished with this kit demonstrates a 2⁺ to 4⁺ intensity reaction. The vial has a listed titer to use as an additional check for the test system (see 1⁺ Dilution Notice). Use this as the calibrator for a 2⁺ to 4⁺ intensity reaction on your microscope.
- Use the HHV6 antibody negative control furnished with this kit as the calibrator for a negative reaction on your equipment.
- Each day's test should include one PBS well in place of a test sample. This is a conjugate control to ensure the conjugate is not reacting with the cell substrate.

1⁺ Dilution Notice

The positive control in this kit is ready to use and provides a 2⁺ to 4⁺ fluorescent intensity when tested. To obtain a 1⁺ fluorescent intensity make two-fold dilutions to the titer indicated on the vial included in this kit. Titer the positive control with the initial use of the kit. The titer you obtain in testing may differ from the listed endpoint titer due to a number of technical reasons. It is best to test the titer indicated on the vial, as well as the two-fold dilution immediately preceding and following the listed titer. It is normal for results obtained for an endpoint (1⁺) titer to differ between laboratories due to factors affecting the intensity of fluorescence. These factors include:

- the power rating of the UV light source in the microscope
- the kind of light source
- the age of the lamp
- the length of the optical path of the microscope and the types of optical filters used
- the accuracy of dilution techniques and the dilution equipment

Limitations of the Procedure

- A serological test such as the IFA serves as an aid to detect viral infection, but its use should not be the sole criteria. The test results should be used in conjunction with information available from the patient, clinical evaluation and other available diagnostic procedures.
- Nonspecific positive reactions such as antinuclear antibody and/or anticytoplasmic antibody reactions can occur in samples from patients with certain autoimmune diseases. Both infected and uninfected cells will fluoresce, and this may obscure a positive HHV6 reaction. Therefore, observation of an autoimmune reaction cannot eliminate the possibility of HHV6 infection.
- SCIMEDX recommends pretreatment of test samples to remove IgG antibody. This additional step helps eliminate false negative and false positive results. When IgG antibody competes with IgM antibody for specific binding sites, IgG antibody can cause a false negative result. When IgG antibody forms immune complexes with the antigenic substrate that may then bind rheumatoid factor (IgM class), IgG antibody can cause a false positive result.
- A sample obtained too early during infection may not contain detectable levels of IgM antibody. If a viral infection is suspected, a second sample should be obtained 7-14 days later. The second sample should be tested concurrently with the first specimen to look for seroconversion or a significant rise in titer of viral-specific IgM or IgG. Either seroconversion or significant rise in titer is indicative of primary infection.
- Because of the possibility of contamination of cord blood with maternal IgM, it is prudent to confirm positive viral IgM antibody results on cord blood samples by testing a follow-up specimen from the infant, preferably within the first five days of life (14).
- Specific IgM antibodies are usually detected in patients with recent primary infection. IgM antibodies may be found in patients with reactivated or secondary infections, or in patients with no other detectable evidence of recent infection (15).

Performance Characteristics

Relative Sensitivity and Specificity: The SCIMEDX HHV6 IgM IFA kit was evaluated relative to a commercially available HHV6 IgM ELISA. The study tested nineteen frozen samples at a commercial clinical laboratory. The overall agreement was 18/19 or 94.7%. The table below summarizes the data.

SCIMEDX IFA				
Alternate ELISA	HHV6 IgM	Reactive	Non-reactive	Total
	Reactive	5	0	5
	Non-reactive	1	13	14
	Total	6	13	19

Reproducibility: Four seropositive sera with various titers (1:10 to 1:80) and six seronegative sera were serially two-fold diluted and assayed three times each on three different assays and the endpoint determined. All thirty endpoint titers were within specifications of $\pm$ one two-fold dilution.

Identical Titer 28/30
 $\pm$ one, two-fold dilution 2/30

Specificity: Thirteen sera with antibody activity to diseases that have the potential for cross-reactivity to HHV6 IgM were tested on this IFA kit. No cross-reactivity was observed for any of the samples. Refer to the table below.

Cross Reactivity Data-SCIMEDX HHV6 IgM IFA

Disease Type	Total Samples	Positive Result
Cytomegalovirus	4	0/4
Epstein-Barr virus	3	0/3
Herpes Simplex Virus 1	4	0/4
Varicella Zoster Virus	2	0/2
Total	13	0/13

Real Time Stability: Real time stability testing of kit components was carried out at 6 month intervals for a minimum of 24 months. The endpoint titer of the positive and negative controls was compared to the endpoint titers at release. Acceptance criteria are endpoint titers within one two-fold dilution of each other. These results were within specifications.

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